

Lakeland Laboratory Supply Order Form

Note: Federal regulations enable us to supply only those materials necessary to submit specimens to our facility for testing.
To request supplies, fax completed supply order form to (269) 983-1904

Date: _____ Client Name: _____

Client Address: _____

Qty Request	Qty Sent	Phlebotomy Supplies	Qty Request	Qty Sent	Cytology/Histology Supplies
		Gold Serum Separator (SST) <input type="checkbox"/> 3.5 <input type="checkbox"/> 5.0			Formalin Biopsy Vial <input type="checkbox"/> 8ml <input type="checkbox"/> 15ml <input type="checkbox"/> 30ml
		<input type="checkbox"/> Green (Lithium Heparin) <input type="checkbox"/> Mint PST			Thin Prep Pap Vials (25 per tray)
		Green (Sodium Heparin)			Thin Prep Brushes/Spatulas (25 each per bag)
		Red Top (serum) Tubes			Thin Prep Brooms (25 per bag)
		Lavender (EDTA) <input type="checkbox"/> 3.0 ml <input type="checkbox"/> 6 ml			CytoLyt® Solution-Fluid Transport Media
		Light Blue (Sodium Citrate)			Microbiology Supplies
		Dark Blue (EDTA)			Anaerobic Transport Tubes
		ACD Solution B Tubes			Routine Culture Swabs
		Pearl Top Tubes			Nasopharyngeal Culture Swabs
		Vacutainer Needles 21G, 1 inch (48 per box)			Stool Culture Collection Kit - Instructions
		Vacutainer Tube Holders (250 per bag)			Ova and Parasite Collection Kit - Instructions
		Needle Disposal Sharps Container <input type="checkbox"/> Small (1 quart) <input type="checkbox"/> Large			Chlamydia/GC PCR Collection Kits-M4RT <input type="checkbox"/> Male <input type="checkbox"/> Female
		Blood Culture Vials <input type="checkbox"/> Adult <input type="checkbox"/> Pediatric			GC Culture Transport Media <input type="checkbox"/> Male <input type="checkbox"/> Female
		Biohazard Bags <input type="checkbox"/> Small <input type="checkbox"/> Large			Group B Strep by PCR Swab
		Lakeland Printer Cartridges			Plain Specimen Container, Sterile
		Kit 6000 - 69G8256			Collection Hats
		M410-412 - 4K00199			M4 Viral Transport Media
		HP 1200 - C7115X			RSV/Viral Screens (double swab)
		HP 1300 - Q2613X			Viral Cultures, Chlamydia Cultures
		HP 2015 - Q7553X			Urinalysis Supplies
		Forms			Urine Cup with integrated sampling device
		Laboratory Test Requisition			Urine Transfer Tube (yellow top)
		Pathology/Cytology Test Requisition			Urine Culture Preservative Tube (gray top)
		Supply Order Forms			Urine Transfer Tube (tiger top preservative)
		Tracking Logs			Kit: Straw, 1 gray/culture
		ABNs			Antiseptic Towelette
		Instruction Sheets: _____			24 Hour Urine Containers - Labels -Instructions

Order Phoned in by: _____ Completed Date/Initials: _____

*Asterisked item not available, lab will send when in stock.